

B R A U N

Chair Topper 2002 -wheelchair roof box

55 kg



Mobility Conversions Ltd, Langley Road South, Salford, Manchester, M6 6SY



Tel: 0161-745-9737 or Fax: 0161-743-0190



Email: info@mobilityconversions.co.uk

Web Site: www.mobilityconversions.co.uk


auto ADAPT
www.autoadapt.com

Chair Topper 2002, wheelchair roof box




Stows your cross frame wheelchair automatically into the roof box

- Weighing only 55 kg the Chair Topper wheelchair roof box is suitable for most car models. Including wheelchair the total roof weight amounts to less than 75 kg.
- Connect the folded cross wheelchair to the hoist device. The stowing process in and out of the roof box is automatic and is operated by means of a hand control.
- With Chair Topper your wheelchair will not stain inside the car or take up space from the luggage area.
- Hand control operated. Stows in 30 sec.
- Alluminium frame. Base/top in ABS plastic.
- Low profile override access cover.
- Available in left or right side loading.
- Low friction side tracks.
- Sealed motors.
- No exposed rivets on the exterior.
- Fully enclosed to be weather tight.
- Integrated manual backup system.

55 kg

Maximum wheelchair dimensions:



A: Width Maximum "folded"	B: Height Maximum "apex"	C: Height Maximum "handle"
22.9 cm	87.6 cm	101.6 cm
25.4 cm	86.8 cm	99.9 cm
27.9 cm	85.9 cm	98.3 cm
30.5 cm	84.9 cm	96.7 cm
33 cm	84.1 cm	95 cm
35.6 cm	83.2 cm	93.3 cm

Maximum wheelchair weight: 20.5 kg

Mobility Conversions Ltd, Langley Road South, Salford, Manchester, M6 6SY



Tel: 0161-745-9737 or Fax: 0161-743-0190

Email: info@mobilityconversions.co.ukWeb Site: www.mobilityconversions.co.uk